HALO2 - Wireless pH Testers



Turn your phone or tablet into a laboratory grade pH Meter.

Bottom line

Hanna's HALO2 Wireless pH Testers combine the best of both worlds, a pH electrode paired with Bluetooth®5.0 Wireless Technology. Take laboratory grade pH measurements using a compatible Apple or Android smart device and the Hanna Lab App. Simple and convenient, test literally anywhere! These testers are perfect for collecting measurements out in the field, in a laboratory or classroom environment, or even as part of the food production process. With 15 styles to choose from for various industries, matching a HALO2 Wireless pH Tester to your testing needs is easy.



Benefits

A Top-of-the-Line Testing Solution

HALO2 Wireless pH Testers are equipped with Hanna's professional line of pH electrodes for industry specific applications. Get the versatility you want and need at a reasonable price! With the HALO2 Line, testing with a professional electrode is no longer limited to just a laboratory setting. Its compact design and data logging capabilities enable you to achieve high quality results from wherever you are.

Calibration Capabilities

Easily perform a 3-point calibration **without** being connected to the Hanna Lab App and up to a 5-point calibration **with** the Hanna Lab app. In addition, set alarm thresholds and calibration reminders via the app.

Trace, Sort, and Share Data

Log pH and temperature measurements every second, tag and annotate measurements, and export into a PDF or CSV format. Visualize and view your measurements on the Hanna Lab App in a graph or table format. Ideal for those performing manual titrations or pH adjustments, all data is tagged with Good Laboratory Practice (GLP) information. This includes calibration buffers used, slope, offset, and date.

Optional Open API Integrations

Optional open API integrations for third-party software are available to reduce data entry. The integration potential is endless whether you want to use it for task management software, laboratory notebooks, etc. Just ask us about it!



When used with HALO2, the Hanna Lab App turns a compatible smart device into a full-featured pH meter. Functions include calibration, measurement, graphing, datalogging, and data sharing as well as Good Laboratory Practice (GLP) capabilities. Measurement and logging start at one second intervals as soon as HALO2 is connected.



HALO2 - Wireless pH Testers



A pH tester for every application. Find the one that works best for you!



HI9810412 Wireless pH Tester for Lab with general purpose electrode



HI9810402
Wireless Refillable pH
Tester for Lab
with general purpose electrode



HI9810422 Wireless pH Tester for Field with specialized electrode



HI9810432 Wireless Refillable pH Tester for Cosmetic Creams with specialized electrode



HI9810302 Wireless pH Tester for Soil with specialized electrode



HI9810332 Wireless Refillable pH Tester for Wine with specialized electrode



HI9810312 Wireless pH Tester for Beer with specialized electrode



HI9810362 Wireless pH Tester for Meat with specialized electrode



HI9810392 Wireless pH Tester for Chocolate with specialized electrode



HI9810382 Wireless pH Tester for Bread & Dough with specialized electrode



HI9810322 Wireless pH Tester for Cheese with specialized electrode



HI9810342 Wireless pH Tester for Milk with specialized electrode



HI9810352 Wireless pH Tester for Sushi with specialized electrode



HI9810372 Wireless pH Tester for Skin & Scalp with specialized electrode



HI9810442 Wireless pH Tester for Leather & paper with specialized electrode